


Al technologies and their application in health care settings

Prof. Farah Magrabi

IMDRF/GMTA Joint Workshop

12 September 2022



Digital health, the promise of Al



"Medical artificial intelligence is primarily concerned with the construction of Al programs that perform diagnosis and make therapy recommendations."

Clancy and Shortliffe, 1984

Now



Supporting decisions:

- diagnosis
- prescribing
- alerts and reminders
- therapy critiquing and planning
- information retrieval
- image recognition and interpretation

Automation in driving

Copyright @ 2021 SAE International. The summary table may be freely copied and distributed AS-IS provided that SAE International is acknowledged as the source of the content. SAE SAE SAE SAE SAE SAE LEVEL 4" LEVEL 3™ LEVEL 5™ LEVEL O" LEVEL 1" LEVEL 2™ You are driving whenever these driver support features You are not driving when these automated driving are engaged - even if your feet are off the pedals and features are engaged - even if you are seated in What does the "the driver's seat" you are not steering human in the driver's seat When the feature You must constantly supervise these support features; These automated driving features have to do? you must steer, brake or accelerate as needed to requests. will not require you to take maintain safety over driving you must drive

machine

human

These are driver support features

These are automated driving features

What do these features do?

These features are limited to providing warnings and momentary assistance These features provide steering OR brake/ acceleration support to the driver

These features provide steering AND brake/ acceleration support to the driver

These features can drive the vehicle under limited conditions and will not operate unless all required conditions are met

This feature can drive the vehicle under all conditions

Example Features

 automatic emergency braking

 blind spot warning

 lane departure warning •lane centering OR

 adaptive cruise control lane centering
 AND

 adaptive cruise control at the same time

traffic jam

chauffeur

taxi
• pedals/
steering
wheel may or
may not be
installed

· local driverless


 same as level 4, but feature can drive everywhere in all conditions






Human information processing





Automation of human information processing



Automating stages of human information processing (Parasuraman et al. 2000)

Open access Original research

BMJ Health & Care Informatics

How machine learning is embedded to support clinician decision making: an analysis of FDA-approved medical devices

David Lyell , Enrico Coiera, Jessica Chen, Parina Shah, Farah Magrabi

To cite: Lyell D, Coiera E, Chen J, et al. How machine learning is embedded to support clinician decision making: an analysis of FDAapproved medical devices. BMJ Health Care Inform 2021;28:e100301. doi:10.1136/ bmihci-2020-100301

Additional material is published online only. To view please visit the journal online (http://dx.doi.org/10.1136/ bmjhci-2020-100301).

Received 07 December 2020 Revised 10 February 2021 Accepted 02 March 2021

ABSTRACT

Objective To examine how and to what extent medical devices using machine learning (ML) support clinician decision making.

Methods We searched for medical devices that were (1) approved by the US Food and Drug Administration (FDA) up till February 2020; (2) intended for use by clinicians; (3) in clinical tasks or decisions and (4) used ML. Descriptive information about the clinical task, device task, device input and output, and ML method were extracted. The stage of human information processing automated by MLbased devices and level of autonomy were assessed. Results Of 137 candidates, 59 FDA approvals for 49 unique devices were included. Most approvals (n=51) were since 2018. Devices commonly assisted with diagnostic (n=35) and triage (n=10) tasks. Twenty-three devices were assistive, providing decision support but left clinicians to make important decisions including diagnosis. Twelve automated the provision of information (autonomous information), such as quantification of heart ejection fraction, while 14 automatically provided task decisions like triaging the reading of scans according to suspected

Summary


What is already known?

- Machine learning (ML)-based clinical decision support (CDS) operates within a human-technology system.
- Clinician interaction with CDS influences how they make decisions affecting care delivery and patient safety.
- Little is known about how emerging ML-based CDS supports clinician decision making.

What does this paper add?

- ML-based CDS approved by the FDA typically provide clinicians with decisions or information to support their decision making.
- Most demonstrate limited autonomy, requiring clinicians to confirm information provided by CDS and to be responsible for decisions.
- We demonstrate methods to examine how MLbased CDS are used by clinicians in the real world.










- Triage (n=10)
- Medical procedures (n=2)
- Treatment (n=1)
- Monitoring (n=1)



- Diagnosis (n=35)
- Triage (n=10)
- Medical procedures (n=2)
- Treatment (n=1)
- Monitoring (n=1)

MACQUARIE University


AUSTRALIAN INSTITUTE OF HEALTH INNOVATION

Image data (n=42)

- CT (n=15)
- MRI (n=10)
- X-ray (n=5)
- Digital breast tomosynthesis (n=3)
- Digital mammography (n=3)
- Echocardiography (n=3)
- Fluoroscopy (n=1)
- Fundus imaging (n=1)
- OCT (n=1)
- PET (n=1)
- Ultrasound (n=1)

Signal data (n=7)

- EKG (n=3)
- Phonocardiography (n=2)
- Polysomnography (n=1)
- Blood glucose & insulin pump (n=1)
- Biometric data (n=1)



- Diagnosis (n=35)
- Triage (n=10)
- Medical procedures (n=2)
- Treatment (n=1)
- Monitoring (n=1)




AUSTRALIAN INSTITUTE OF HEALTH INNOVATION

Image data (n=42)

- CT (n=15)
- MRI (n=10)
- X-ray (n=5)
- Digital breast tomosynthesis (n=3)
- Digital mammography (n=3)
- Echocardiography (n=3)
- Fluoroscopy (n=1)
- Fundus imaging (n=1)
- OCT (n=1)
- PET (n=1)
- Ultrasound (n=1)

Signal data (n=7)

- EKG (n=3)
- Phonocardiography (n=2)
- Polysomnography (n=1)
- Blood glucose & insulin pump (n=1)
- Biometric data (n=1)



- Quantification (n=13)
- Triage notification (n=10)
- Case level finding of disease (n=6)
- Identify features of disease (n=6)
- Clinical grades (n=5)
- Enhanced images (n=4)
- Automatic coding of features (n=2)
- Automatic control of devices (n=2)
- Treatment recommendations (n=1)

- Diagnosis (n=35)
- Triage (n=10)
- Medical procedures (n=2)
- Treatment (n=1)
- Monitoring (n=1)




AUSTRALIAN INSTITUTE OF HEALTH INNOVATION

Image data (n=42)

- CT (n=15)
- MRI (n=10)
- X-ray (n=5)
- Digital breast tomosynthesis (n=3)
- Digital mammography (n=3)
- Echocardiography (n=3)
- Fluoroscopy (n=1)
- Fundus imaging (n=1)
- OCT (n=1)
- PET (n=1)
- Ultrasound (n=1)


Signal data (n=7)

- EKG (n=3)
- Phonocardiography (n=2)
- Polysomnography (n=1)
- Blood glucose & insulin pump (n=1)
- Biometric data (n=1)



- Artificial intelligence (n=15)
- Machine learning (n=14)
- Deep learning (n=11)
- Convolutional neural network (n=6)
- Neural network (n=1)
- Deep neural network (n=1)
- Deep convolutional neural network (n=1)


- Quantification (n=13)
- Triage notification (n=10)
- Case level finding of disease (n=6)
- Identify features of disease (n=6)
- Clinical grades (n=5)
- Enhanced images (n=4)
- Automatic coding of features (n=2)
- Automatic control of devices (n=2)
- Treatment recommendations (n=1)






AUSTRALIAN INSTITUTE OF HEALTH INNOVATION


Level of autonomy



1. Al provides information for humans to use in decision-making




Al-Rad (Siemens Medical Solutions)



Quantification of

- lung volume
- lung lobes
- specified lung lesions

Autonomous Information




2. Al assists human decisions



Cardiologs (Cardiologs Technologies)


"Clinicians must not use the [AI] generated output as the primary interpretation."



Identification of abnormal cardiac rhythms

responsibility for final decision resides with clinician

Assistive




3. Al decides in place of human expert



IDx-DR (Idx LLC)

Autonomous Decision






n = 14

Safely integrating Al into clinical workflow










Summary



- Framework to classify AI/ML-based clinical decision support by level of autonomy
- Current systems:
 - provide information
 - <u>assist humans</u> to make decisions
 - provide <u>decisions</u>
- Last mile of safe implementation and use needs careful attention






AUSTRALIAN INSTITUTE OF HEALTH INNOVATION


Thank you

farah.magrabi@mq.edu.au

Assistive




Autonomous Information




Autonomous Decision





Automation of human information processing



Automating stages of human information processing (Parasuraman et al. 2000)